

TERMINATION:-CLAMPING SLEEVES:-WARNING: NOMINAL OPERATING VOLTAGE 210KV DC STANDARD FLANGE T20172/A CABLE SHOULD ONLY BE FITTED BY A TRAINED PERSON.
OIL MUST NOT BE USED ON THE RUBBER TERMINATION, ONLY SUITABLE HIGH VOLTAGE GREASE. - USING AN ESSEX X-RAY APPROVED RECEPTACLE CABLE LENGTH TOLERANCE +/- 100mm QUICK LOCK T20173 ENSURE BY ADJUSTING THE CLAMPING SLEEVE THAT WHEN THE TERMINATION IS FULLY INSERTED IN THE SOCKET THERE IS A CABLE:- VOLTAGE WITHSTAND GAP OF 6mm BETWEEN THE MATING FACES OF THE CLAMPING SLEEVE 50KV (RMS) AC INSTALLATION GREASE P00075 AND SOCKET BEFORE TIGHTENING THE CLAMPING SLEEVE. FOR FULL CABLE DATA SEE ESSEX CABLE DATA SHEET A

REVISIONS

ESSEX

Dandong Light Electric Co., Ltd Rm.1606 Jiadi Square, No.64, Binjiang Middle Rd, Zhenxing, Dandong, Liaoning China T:0415-6123232 F:0415-6123533

NO.	DATE	BY	DRAWN	APPROVED	DATE	ISSUE	SCALE		© NOT TO BE REPRODUCED IN WHOLE OR PART WITHOUT
6	11.05.01	DCN 1455	JDF		15.09.95	9		N.T.S.	PRIOR PERMISSION OF ESSEX X-RAY & MEDICAL EQUIPMENT LTD.
7	13.07.01	DCN 1534	TITLE	, 05,55		IT AND	. –	(D70	DRAWING No.
8	08.07.04	DCN 2047	210K\	/ SEIFEF	RT RIGH	HT ANG	JYDE	(R30	INCIONE
9	06.03.09	DCN 2905]		((PARLE	IYPE	C20	42)

DUE TO CONTINUING PRODUCT DEVELOPMENT ESSEX X-RAY RESERVES THE RIGHT TO MODIFY DATA WITHOUT NOTICE