

THIS DATA SHEET IS FOR IDENTIFICATION PURPOSES ONLY. NOMINAL DIMENSIONS IN MM.

TECHNICAL DATA:	ACCESSORIES:	INSTALLATION:
TERMINATION:- NOMINAL OPERATING VOLTAGE 65KV DC - USING AN ESSEX X-RAY APPROVED RECEPTACLE	CLAMPING SLEEVES:- FLANGE T20051/A	WARNING: CABLE SHOULD ONLY BE FITTED BY A TRAINED PERSON.
CABLE LENGTH TOLERANCE +/- 100mm CABLE:- VOLTAGE WITHSTAND 100KV DC 30KV (RMS) AC 🚖	FLANGE T20097 INSTALLATION GREASE P00075	OIL MUST NOT BE USED ON THE RUBBER TERMINATION, ONLY SUITABLE HIGH VOLTAGE GREASE. ENSURE BY ADJUSTING THE CLAMPING SLEEVE THAT WHEN THE TERMINATION IS FULLY INSERTED IN THE SOCKET THERE IS A GAP OF 4mm BETWEEN THE MATING FACES OF THE CLAMPING SLEEVE AND SOCKET BEFORE TIGHTENING THE CLAMPING SLEEVE.
FOR FULL CABLE DATA SEE ESSEX CABLE DATA SHEET	INSTALLSTION SKENE 1999/9	AND SUCKET BEFORE HIGHTEINING THE CLAMPING SLEEVE.

ESSEX

Dandong Light Electric Co.,Ltd Rm.1606 Jiadi Square,No.64,Binjiang Middle Rd, Zhenxing, Dandong, Liaoning China T:0415-6123232 F:0415-6123533

	REVISION:	S	DUE TO C	DUE TO CONTINUING PRODUCT DEVELOPMENT ESSEX X-RAY RESERVES THE RIGHT TO MODIFY						DATA WITHOUT NOTICE	
NO.	DATE	BY	DRAWN	APPROVED	DATE	ISSUE	SCALE		O NOT TO BE REPRODUCED	IN WHOLE	OR PART WITHOUT
3	02.01.01	DCN 1405	JDF		07.08.96	6		N.T.S.	PRIOR PERMISSION OF ES	SSEX X-RAY	/ & MEDICAL EQUIPMENT LTD.
4	11.05.01	1400	TITLE 65KV TERMINATION (CABLE TYPE C2212)							DRAWING No.	
5	05.06.06	DCN 2482								DS1010	
6	06.03.09	DCN 2905	(CADLL TIPL CZZTZ)								טוטוטן