

TECHNICAL DATA:	ACCESSORIES:	INSTALLATION:
TERMINATION:— NOMINAL OPERATING VOLTAGE 225KV DC - USING AN ESSEX X—RAY APPROVED RECEPTACLE CABLE LENGTH TOLERANCE +/- 100mm CABLE:— VOLTAGE WITHSTAND 250KV DC	CLAMPING SLEEVES:— STANDARD FLANGE T20289 INSTALLATION GREASE P00075	WARNING: CABLE SHOULD ONLY BE FITTED BY A TRAINED PERSON. OIL MUST NOT BE USED ON THE RUBBER TERMINATION, ONLY SUITABLE HIGH VOLTAGE GREASE. ENSURE BY ADJUSTING THE CLAMPING SLEEVE THAT WHEN THE TERMINATION IS FULLY INSERTED IN THE SOCKET THERE IS A GAP OF 6mm BETWEEN THE MATING FACES OF THE CLAMPING SLEEVE
75KV (RMS) AC 🛕 FOR FULL CABLE DATA SEE ESSEX CABLE DATA SHEET 🛕		AND SOCKET BEFORE TIGHTENING THE CLAMPING SLEEVE.

ESSEX

Dandong Light Electric Co.,Ltd Rm.1606 Jiadi Square,No.64,Binjiang Middle Rd, Zhenxing, Dandong, Liaoning China T:0415-6123232 F:0415-6123533

REVISIONS		DUE TO C	ONTINUING P	RODUCT DE	EVELOPMEN	T ESSEX	X-RAY	RESERVES	THE RIGI	OT TH	MODIFY	DATA WI	THOUT NOTI	CE	
NO.	DATE	BY	DRAWN	APPROVED	DATE	ISSUE	SCALE		С мот то	BE REPRO	DUCED	IN WHOLE	OR PART V	VITHOUT	
3	02.01.01	DCN 1405	JDF		26.03.97	6		N.T.S.	PRIOR I	PERMISSION	OF ESS	SEX X-RA	/ & MEDICA	L EQUIPMENT	LTD.
4	11.05.01	DCN 1455	TITLE	TILE 225KV ESSEX STRAIGHT (R28)								DRAWING	No.		
5	13.06.02	DCN 1701												100	Λ
6	06.03.09	DCN 2905	(CABLE TYPE C2236)									102	+		